

Ministry of Labour, Employment & Social Security

National Office of Equipments & Accessories for Persons with Disabilities



ONAAPH

PRODUCT CATALOGUE

WALKING AIDS

**For the full inclusion of persons with special needs
in social and professional life**



www.onaaph.dz



WALKING AIDS



ONAAPH offers a wide range of technical mobility aids designed to enhance movement and support individuals experiencing difficulties with ambulation. The following is an overview of available products:



Electric Wheelchair

Indicated for patients with conditions that prevent them from walking or self-propelling a manual wheelchair. Conditions may be of a medical nature (myopathy, poliomyelitis...) or of a traumatic nature (spinal cord injuries).



ALUMINIUM WHEELCHAIR

This aluminum wheelchair is lightweight and robust, indicated for patients who are able to self-propel the wheelchair.

REINFORCED STEEL -WHEELCHAIR

Designed for users weighing over 90 kg. Combines robustness and comfort to meet the specific needs of persons requiring reinforced support.



COMMODE WHEELCHAIR

Consists of a welded tube frame designed to accommodate a body support system and a toilet seat (WC seat), providing hygienic comfort and practicality for daily care.



ACTIVE WHEELCHAIR

Designed to meet the needs of active and demanding users who travel extensively on their own or engage in sustained activity, it is equipped with a low backrest to allow a great range of motion for the upper limbs.

Wheelchairs for Cerebral Palsy (IMC)

IMC WHEELCHAIR



Designed to provide adapted postural support, with multiple adjustments for comfort and stability. It enables facilitated mobility while ensuring safety and ergonomics to meet the specific needs of users.

IMC STROLLER (Poussette IMC)

Designed for children with motor disorders, offering comfort, postural support and safety. It is adjustable and tailored to the specific needs of children with cerebral palsy.



CANES & CRUTCHES



T-CANE / C-CANE

C-canes have a curved handle for light support.

T-canes offer a better grip for increased stability.



TRIPOD CANE


Offers better stability thanks to its three-point base. Ideal for persons with balance problems.



CANADIAN CRUTCH (Forearm Crutch)


Provides a forearm support for better weight distribution and increased stability.

Ideal for persons in rehabilitation or with mobility difficulties.



CRUTCH (Axillary / Underarm)

Relieves an injured leg by providing support to facilitate walking and maintain balance.



WHITE CANE

Designed for visually impaired or blind persons, allowing them to detect obstacles and move around safely and autonomously.

ROLLATOR

A four-wheeled walking aid with brakes and a built-in seat, offering stability and support while allowing the user to rest when needed.



WALKER (D»ambulateur)

Equipped with wheels on the front feet, facilitating progression and reducing the effort required to move. Provides stable support while allowing better maneuverability.



WALKING FRAME (Cadre de Marche)

A four-legged support device that folds for easy storage and transport. Offers additional stability and support while being practical to carry.





PRODUCT CATALOGUE

**National Office of Equipements and Accessories
for Persons with Disabilities**



City Ben Boulaid, Bourouba, Algiers



SCAN QR

www.onaaph.dz

